MS303-1 Rev 3
11-2019

Labeling Symbols:
 Symbols may be used on some international package labeling for easy identification.

Caution, consult accompanying documents

Catalog Number

USA law restricts this device to sale by or on the order of physician

Manufacturer

Expiration Date

Authorized European Representative

Do not reuse after opening

Do not use if the product sterilization barrier or its packaging is compromised.

Lot number

Method of sterilization – gamma

Temperature Limitation

EC Representative
 MDSS Medical Device Safety Service GmbH
 Schiffgraben 41
 30175 Hannover, Germany

Manufacturer:
 Collagen Matrix, Inc.
 15 Thornton Road
 Oakland, NJ 07436 USA

Zimmer® Socket Repair Membrane
 Instructions For Use

Distributed by:
 Zimmer Biomet Dental
 4555 Riverside Drive
 Palm Beach Gardens
 FL 33410 USA
 +1.561.776.6700
 zimmerbiomet.com

Intended Use: Zimmer® Socket Repair Membrane is a bioresorbable, implantable collagen material that is intended for use in dental surgery procedures as a material for placement in the area of dental implant, bone defect or ridge reconstruction to aid in wound healing post dental surgery.

Description: Zimmer® Socket Repair Membrane is a white, nonfriable membrane matrix engineered from highly purified type I collagen fibers derived from bovine Achilles tendon. Zimmer® Socket Repair Membrane is resorbable which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane.

Zimmer® Socket Repair Membrane has a morphology of dense oriented fibers for mechanical strength. Macromolecular permeation studies have shown that the membrane is permeable to macromolecules. Its porosity is such that it effectively retards epithelial down growth and prevents gingiva connective cell migration into the wound site. The semi-permeability properties of the membrane permit the exchange of essential nutrients for wound healing.

Zimmer® Socket Repair Membrane is supplied sterile, non-pyrogenic and for single use only.

Administration: Zimmer® Socket Repair Membrane is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The periodontal defect or bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). The clinician should perform thorough debridement and good planning of the defect. Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect. As much tissue as possible should be preserved to allow for primary closure of the wound and correct positioning of the flaps.

Zimmer® Socket Repair Membrane can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to the final placement.

Zimmer® Socket Repair Membrane can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

Zimmer® Socket Repair Membrane should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. An attempt should be made to suture the mucoperiosteal flap over the collagen membrane to avoid accelerated resorption due to membrane exposure.

Post-operative Procedures:

Zimmer® Socket Repair Membrane is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior for four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

Zimmer® Socket Repair Membrane should be completely resorbed 26 to 38 weeks following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Contraindications:

Zimmer® Socket Repair Membrane is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity
- known allergy to collagen of animal origin or other bovine-derived products
- clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing.

Warning:

Clinicians should use care in screening their patients for any known allergies to collagen or bovine-derived products. Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to Zimmer® Socket Repair Membrane.

Precautions:

As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, including a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little periodontium or bone.

Zimmer® Socket Repair Membrane cannot be re-sterilized. Open, unused Zimmer® Socket Repair Membrane must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

Adverse Reactions:

Possible complications that can occur with any dental surgery include infection, swelling of the introral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the use of anesthesia. Minor discomfort may occur for a few days.

Safety:

The product is manufactured from bovine Achilles tendon, which is classified as tissues with no detected infectivity for Bovine Spongiform Encephalopathy, BSE (World Health Organization Guidelines). The bovine tendon is known to be one of the richest sources of type I collagen that is commercially available.

The manufacturing process for the product meets European Standards and International Standards for animal tissue sourcing, handling and inactivation of Spongiform Encephalopathy (SE) pathogens. This process involves a treatment with sodium hydroxide that is a recognized method of inactivation of SE pathogens.

A viral inactivation study for the product's manufacturing process was conducted by an independent laboratory. In this study, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the following viral strains: Bovine Viral Diarrhea (enveloped virus) and Porcine Parvoviridae (non-enveloped virus). The study results showed that each of the manufacturing steps evaluated, including the sodium hydroxide treatment, is effective in inactivating these viruses.

Storage:

The product should be stored at room temperature. Avoid excessive heat and humidity.

How Supplied:

One (1) membrane per package, size 10mm x 20mm

Caution:

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

MS303-1 Rev 3
11-2019

Caution, consult accompanying documents

Catalog Number

USA law restricts this device to sale by or on the order of physician

Manufacturer

Expiration Date

Authorized European Representative

Do not reuse after opening

Do not use if the product sterilization barrier or its packaging is compromised.

Lot number

Method of sterilization – gamma

Temperature Limitation

EC Representative
 MDSS Medical Device Safety Service GmbH
 Schiffgraben 41
 30175 Hannover, Germany

Manufacturer:
 Collagen Matrix, Inc.
 15 Thornton Road
 Oakland, NJ 07436 USA

Surgical Procedure and Technique
Step 1 – Atraumatically Extract Tooth Using a Flapless Technique: Prepare the bone graft material (e.g. Puros® Particulate) according to the instructions for use. Utilizing an atraumatic, flapless technique, carefully extract the tooth and completely debride the socket.

Step 2 – Assess Shape & Prepare Membrane: Remove the Zimmer® Socket Repair Membrane from the packaging. Insert the small end of the pre-cut membrane into the socket, ensuring that it extends over the facial defect both laterally and apically. Establish how much of the membrane's large end will be needed to fold closely over the top of the socket to facilitate closure after grafting. Remove the membrane and trim as needed.

Step 3 – Line Facial Defect with Membrane: Insert the small end of the membrane into the socket and against the facial tissue.

Caution: To help maintain blood supply to the remaining facial plate, do not lift the periosteum off the bone.

Step 4 – Augment Socket with Bone Graft Material: Augment the socket with graft material using a sterile syringe or sterile applicator. To fill the syringe, pull back on the plunger and gently press its tip into the hydrated graft material. Repeat until the plunger is filled. After placing bone graft material into socket, firmly compress with a sterile instrument. This will prevent voids at the apex of the socket and push the facial tissue labially for better ridge shape.

Step 5 – Close Socket with Membrane: Gently close the socket by folding the membrane over the bone graft material until the membrane approximates with the palatal/lingual tissue. The membrane should be stable due to the graft material pressing it against the facial plate.

Step 6 – Suture Membrane to Soft Tissue: Place at least 2 to 3 interrupted 5-0 absorbable sutures to secure the membrane to the palatal/lingual tissue. Note: Typically, no sutures on the labial/buccal aspect are needed since the membrane will be held by pressure of the graft material.

Zimmer is a registered trademark of Zimmer Biomet or its affiliates.

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Rappresentante CE: MDS5 Medical Device Safety Service GmbH Schifffgraben 41 30175 Hannover, Germania
Produttore: Collagen Matrix, Inc. 15 Thornton Road Oakland, NJ 07436 USA
PEU-representant: MDS5 Medical Device Safety Service GmbH Schifffgraben 41 30175 Hannover, Tyskland
Fabrikant: Collagen Matrix, Inc. 15 Thornton Road Oakland, NJ 07436 USA