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CopiOs® Cancellous Particulate Xenograft

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Please read carefully!

CopiOs® Cancellous Particulate Xenograft

Regenerative matrix for bone regeneration

STERILE – single use for one patient only

Product name

CopiOs Cancellous Particulate Xenograft

Description

CopiOs Cancellous Particulate Xenograft is a matrix for bone regeneration, made from solvent-preserved bovine cancellous bone tissue sterilised by gamma irradiation.

The tissue is selected according to the strictest specifications. CopiOs Cancellous Particulate Xenograft is preserved using the company's proprietary multiple step process (Tutoplast™ process), which consists of a thorough purification of the tissue and gentle solvent dehydration; the process is subject to continuous quality controls. Preservation also includes sterilisation of the finished product with low-dose gamma irradiation. The Tutoplast™ process maintains the natural collagen structure and thus the high biomechanical stability of the tissue.

Mechanism of action

CopiOs Cancellous Particulate Xenograft functions as a template for the local regeneration of vascularised bone tissue. Absorption and tissue regeneration processes usually start a few days postoperatively and last for several weeks or even months. Bone remodelling depends on the size and volume of the CopiOs bone particulate, and the capacity for regeneration of the implant bed.

Indications for use

CopiOs Cancellous Particulate Xenograft is intended for use in oral surgical procedures as a resorbable material for the following indications:

- Regeneration of periodontal bone defects
- Regeneration of cavities between the alveolar wall and immediate implants
- Horizontal alveolar ridge augmentation
- Alveolar ridge augmentation at implant sites with sufficient residual bone and a good blood supply

Contraindications

- Known hypersensitivity to bovine collagen or bovine trabecular bone
- Disorders or diseases carrying an unacceptable increase in postoperative risk
- Disorders or conditions which inhibit healing due to poor vascularisation, like heart disease, diabetes or smoking
- Active or latent infection in the area of the graft

Adverse reactions

Local reactions including inflammation, granuloma, foreign body reaction or other types of adverse events including slow resorption or infections may occur in patients who are not known to be hypersensitive to foreign collagen. These adverse reactions are typically self-limiting and are generally expected to resolve rapidly when treated appropriately.

Patients should be told to inform the doctor about any possible side effects. All adverse reactions must be reported to Tutogen Medical GmbH.

Warnings

CopiOs Cancellous Particulate Xenograft should be implanted adjacent to well-vascularised tissue for best infiltration of cells, revascularisation and incorporation.

CopiOs Cancellous Particulate Xenograft should not be re-sterilised or re-used and should be discarded if the sterile packaging is opened or damaged.

CopiOs Cancellous Particulate Xenograft remains sterile, as long as the sterile packaging is not opened or damaged. Sterility is no longer ensured once the pack has been opened. If the product is re-sterilised, the stability/quality of the material can no longer be guaranteed. Re-use is therefore not permitted as it may harm the patient's health.

Integration and remodeling can be affected by comorbidities or conditions that lead to a reduced vascularisation, such as diabetes or smoking.

Precautions for use

Limited clinical data suggest CopiOs Cancellous Particulate Xenograft may be used safely in infected, contaminated, or possibly infected or contaminated areas; however, the clinical data are not sufficient to support routine use of CopiOs Cancellous Particulate Xenograft in infected, contaminated or possibly infected or contaminated areas and such uses should be weighed carefully with a benefit-risk assessment for each individual patient.

The Tutoplast™ process is highly effective against all pathogenic species. However, as with all biological products, the transmission of infectious diseases can be ruled out only to the extent of current knowledge.

Inspect the device, including all packaging and labeling materials carefully:

- Do not use past expiration date.
- Do not use if the device or packaging is damaged.
- Do not use if any discrepancies in label information are found.
- Do not use if device was rehydrated in any solution above 37°C.
- Dispose of unused CopiOs Cancellous Particulate Xenograft properly.

Before the clinical use of CopiOs Cancellous Particulate Xenograft, the surgeon should be familiar with all the aspects and limitations of the surgical techniques and grafting. Patients must be informed about the limitations of graft procedures.

Handling instructions

Follow general rules for aseptic handling of medical devices when using CopiOs Cancellous Particulate Xenograft.

Preparation

- Select appropriate CopiOs Cancellous Particulate Xenograft quantity/dimensions for defect to be filled.
- Remove the pouch from the carton and peel open the pouch.
- Remove the vial using aseptic technique.
- Before implantation, rehydrate by covering CopiOs Cancellous Particulate Xenograft with sterile 0.9 % saline solution (maximum temperature 37°C).
- Stir and let stand for at least two minutes.
- Rehydrated CopiOs Cancellous Particulate Xenograft can be kept in this solution throughout surgical procedure until needed.
- Autologous blood, blood components or bone marrow aspirate may be added only after rehydration and immediately before implantation.

Implantation

- Prepare surgical site using standard techniques.
- Insert CopiOs Cancellous Particulate Xenograft firmly into graft site to prevent displacement and allow stable osseointegration of graft material.
- Ensure every effort is made to provide tension-free, hermetic wound closure.
- Properly dispose of unused CopiOs Cancellous Particulate Xenograft.
- Instruct patients to maintain proper dental hygiene as a basic requirement for successful augmentation.

Interactions

None known

Magnetic Resonance (MR) Environment:

CopiOs Cancellous Particulate Xenograft is composed of natural trabecular bone collagen and is safe for the MR environment with no MR artefacts expected.

Shelf life

Do not use CopiOs Cancellous Particulate Xenograft after the expiry date printed on the pouch and outer packaging.

Storage instructions

Store this medical device at room temperature between 15°C and 30°C, in a clean, dry place; protected from direct sunlight.

Note

As a product of natural origin, the colour of CopiOs Cancellous Particulate Xenograft may vary.

Date of last revision

2019/04 Rev04

CopiOs® is a registered trademark of Zimmer Biomet or its affiliates.

Manufacturer:

Tutogen Medical GmbH
 Industriestrasse 6
 91077 Neunkirchen am Brand
 Germany
 Tel.: +49 (0) 9134 99 88-0
 Fax: +49 (0) 9134 99 88-99
 tutogen@rtix.com

Distributed by:

ZIMMER BIOMET
 BIOMET 3i Dental Ibérica, S.L.
 WTC Almeda Park, Ed. 4, Pl. 2
 C/Tirso de Molina, 40
 08940, Cornellà de Llobregat, Spain
 Tel: +34 93 470 59 50
 Fax: +34 93 372 11 25
 zimmerbiometdental.com

	Do not re-use		Manufacturer		Keep away from sunlight
	Do not re-sterilize		Sterilized using irradiation		Keep dry
	Batch code		Use-by date		Temperature limit
	Catalogue number		Consult instructions for use		Do not use if the package is damaged

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