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**Information for use**  
Please read carefully!

## PUROS® Dermis Allograft Tissue Matrix

**STERILE (STERILE R) – for single use only**

**Product name**  
Puros® Dermis Allograft Tissue Matrix

**Description**  
Puros® Dermis is a collagenous membrane derived from solvent preserved,  $\gamma$ -irradiated human Dermis. This tissue is selected in compliance with strict specifications. Puros® Dermis is preserved using a proprietary, multi-step process (Tutoplast® process) which consists of meticulous cleaning of the tissue and gentle solvent dehydration monitored by continuous quality controls. Preservation also includes sterilisation of the finished product by low-dose gamma irradiation. The natural collagen structure and consequently the high biomechanical stability of this tissue are preserved due to the Tutoplast® process.

**Other components**  
None.

**Form and contents**  
Soft Tissue Graft of human origin (Allograft)

**Mode of action**  
For implantation. Puros® Dermis works as a biological scaffold for the new formation of vascularized connective tissue. The process of resorption and new tissue formation generally begins after a few days post surgery and extends over several weeks and months. This remodeling process depends on the size of the tissue and on the remodeling capacity of the recipient site.

**Indications**  
For implantation as a replacement or reinforcement of lost, resected or weakened connective tissue (e.g. Recession Coverage in Regenerative Periodontology or Implantology).

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**Contraindications**

- Active or latent infection in or around the implantation site.
- Any condition or diseases carrying an unacceptable increase in the postoperative risk.

**Handling instructions**

Please note that the general principles of sterile handling should be applied when using Puros® Dermis.

- Open package by firmly grasping both sides of the pouch top and peeling outward. Place the tissue graft onto a sterile field.
- Puros® Dermis should be trimmed dry.
- Afterwards, rehydrate Puros® Dermis with sterile isotonic saline solution or patient's blood until the membrane becomes soft and flexible, allowing the membrane to adapt to the surface very well.
- The graft should be sized according to the defect and securely placed to prevent displacement and aid incorporation. Implant the graft so that outer edges of the Puros® Dermis do not protrude.
- If fixation of Puros® Dermis is selected, either resorbable or non-resorbable suture material (ideally monofilament) may be used with a round trauma needle. Place the stitches 2-3 mm from the edge of the implant where possible.
- A tension-free, hermetic closure of the flap should be attempted.
- In case of dehiscence, the uncovered Puros® Dermis does not have to be explanted categorically. The wound site should be disinfected daily (e.g. with chlorhexidine) until the membrane is covered with mucosa.
- The process of absorption and tissue neogenesis begins after 1-2 days and may last for weeks or months, depending on the size of the graft and reactivity of the host site.

**Interactions with other agents**  
None known.

**Precautions**  
The Tutoplast® process is highly active against all kinds of pathogens and has been used in more than 3.000.000 surgical operations with Tutoplast preserved tissue over the last 35 years. As with all biological products, the transmission of infectious diseases can not be ruled out categorically.

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**Adverse reactions**  
None known.  
Patients should be instructed to inform the physician if adverse reactions should occur. Adverse reactions should be immediately reported to Zimmer Biomet or Tutogen Medical GmbH.

**Expiry date**  
Puros® Dermis has a shelf life of 5 years in the undamaged sterile package. Puros® Dermis must not be used after the expiry date printed on the primary package and on the outer package. If the primary package has been opened, but the Puros® Dermis has not been used, the graft must subsequently be discarded. The tissue graft cannot be re-sterilised.

**Special precautions for storage**  
Store Puros® Dermis in a clean and dry place, protected from direct sun exposure at room temperature between 15°C and 30°C.

**Special precautions for disposing of unused products**  
None

**Date of the package insert**  
06/2017 „02“

**Tutoplast® is a registered trademark of Tutogen Medical GmbH; PUROS® is a registered trademark of Zimmer Biomet or its affiliates.**

Tutoplast® processed allografts have been successfully applied for more than 35 years in over 3.000.000 surgical procedures. This success is based, on one hand, on strict criteria of donor selection and serological tests. On the other hand, Puros® Dermis is preserved with the multi-step Tutoplast® procedure that involves purification of the tissue and solvent dehydration. All manufacturing steps are subject to strict quality requirements and are monitored by continuous quality controls. This process removes cellular components, non-collagenous proteins and, thus, antigenic properties. This process does not result in any chemical residues and is highly active against all types of pathogens. Preservation also includes sterilisation of the finished product by low-

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