# RegenerOss® Resorbable Xenograft Syringe

## Porcine Anorganic Bone Mineral in Syringe

### Instructions for Use

**Manufacturer:**  
Collagen Matrix, Inc.  
15 Thornton Road  
Oakland, New Jersey 07436 USA  
[www.collagenmatrix.com](http://www.collagenmatrix.com)

**Authorized EU Representative:**  
MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany

**Distributed By:**  
Zimmer Biomet Dental  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410 USA  
+1.561.776.6700  
[zimmerbiometdental.com](http://zimmerbiometdental.com)

### ENGLISH

#### RegenerOss® Resorbable Xenograft Syringe

##### Porcine Anorganic Bone Mineral in Syringe

**INSTRUCTIONS FOR USE.**  
A sterile, biocompatible anorganic porous porcine bone mineral for use in periodontal, oral and maxillofacial surgery.

**Description:**  
RegenerOss® Resorbable Xenograft Syringe is a porous bone mineral matrix consisting predominantly of calcium phosphate, preloaded into a delivery applicator. RegenerOss® Resorbable Xenograft is produced by removal of the organic components from porcine cancellous bone. The delivery applicator is an open barrel syringe with a perforated cap for ease of hydration and placement of RegenerOss® Resorbable Xenograft. RegenerOss® Resorbable Xenograft Syringe is sterilized by gamma irradiation. It is supplied sterile, non-pyrogenic and for single use only.

**Properties/Actions:**  
The anorganic bone matrix of RegenerOss® Resorbable Xenograft has interconnecting macro- and microscopic porous structure that supports the formation and ingrowth of new bone at the implantation site. RegenerOss® Resorbable Xenograft may be considered when autogenous bone is not indicated, or insufficient in quantity to fulfill the needs of the proposed surgical procedure.

**Indications and Usage:**  
RegenerOss® Resorbable Xenograft Syringe is indicated for:

- Augmentation or reconstructive treatment of the alveolar ridge.
- Filling of infrabony periodontal defects.
- Filling of defects after root resection, apicoectomy, and cystectomy.
- Filling of extraction sockets to enhance preservation of the alveolar ridge.
- Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR).
- Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR).

**Instructions for Use:**  
After exposure of the bony defect with a mucoperiosteal flap, all granulation tissue must be carefully removed.

- Remove the sterile pre-loaded RegenerOss® Resorbable Xenograft Syringe from packaging and hydrate with osseous coagulum, patient's blood or sterile normal saline, as follows:

- with the perforated cap left on the syringe, immerse the tip of the syringe in the fluid;
- pull up on the plunger to draw the fluid into the syringe;
- push down on the plunger to expel excess fluid/air;
- pull and push action should be repeated three times or until the material is fully hydrated.

- Once hydrated, remove the perforated cap and deliver the hydrated granules to the defect site.

- In order to assure the formation of new bone, RegenerOss® Resorbable Xenograft should only be placed in direct contact with well vascularized bone. Cortical bone should be mechanically perforated.

- If large maxillofacial defects are present, RegenerOss® Resorbable Xenograft should be mixed with autogenous bone in a ratio of approximately 1:1. In this situation, RegenerOss® Resorbable Xenograft should be poured out of the delivery syringe for mixing with autogenous bone and implanted without the use of the delivery syringe applicator. RegenerOss® Resorbable Xenograft is also available in a jar without the delivery syringe.

- Overfilling of the defects should be avoided.

- The mucoperiosteal flaps should be sutured to achieve primary closure, if possible. A surgical dressing may be placed over the wound for one to two weeks.

- If primary wound closure cannot be achieved completely, further immobilization of the flap (e.g., by incision through the periosteum) should be performed and/or a bioabsorbable membrane (e.g. Resorbable Collagen Membrane) should be placed over the bone graft site.

- Post-operative procedures: Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for 2 weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent. Consult the manufacturer's Instructions for Use for side effects and other known risks associated with the use of chlorhexidine gluconate oral rinses.

**Contraindications:**  
Contraindications customary to the use of bone grafts should be observed. RegenerOss® Resorbable Xenograft Syringe should not be used in patients with:

- Acute or chronic infection (osteomyelitis) at the surgical site

- Metabolic diseases (diabetes, hyperparathyroidism, osteomalacia)

- Severe renal dysfunction, severe liver disease

- High dose corticosteroid therapy

- Vascular impairment at the implant site

- Osteoporosis

- Known allergy or hypersensitivity to porcine-derived implantable materials

**Warnings:**  
The device should be secured to prevent motion and migration; use in areas where the graft can be adequately contained.

Do not use if package is opened or damaged or if expiration date has been exceeded.

RegenerOss® Resorbable Xenograft Syringe cannot be re-sterilized or re-used. Open, unused product must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

**Precautions:**  
In order to facilitate the formation of new bone, RegenerOss® Resorbable Xenograft should only be implanted in direct contact with a well vascularized bony tissue. Drilling may be recommended to facilitate bleeding from cortical bone. In larger defects a mixture of autogenous bone or bone marrow may improve the formation of new bone. The use of RegenerOss® Resorbable Xenograft Syringe has not been evaluated in children.

**Implantology:**  
Generally, in augmented areas, the placement of titanium fixtures should take place once the bone has sufficient strength and integrity for dental implant placement, which is typically greater than 6 months

after implantation of a bone graft material. For sinus floor elevation, typically 9-12 months should be allowed after implantation of bone graft material before placement of the titanium fixtures. X-rays should be taken to confirm the bone integrity prior to dental implant placement.

#### Periodontology

The filling of periodontal defects with RegenerOss® Resorbable Xenograft requires (along with plaque control) the successful local treatment of the periodontal lesion (e.g. root planing, debridement of granular tissue) prior to implantation.

#### Caution:

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

#### Adverse Reactions:

Possible complications that may occur with any dental surgery include swelling at the surgical site, flap sloughing, bleeding, local inflammation, bone loss, infection or pain.

#### How Supplied:

RegenerOss® Resorbable Xenograft Syringe is supplied sterile, non-pyrogenic, and for single use only.

Catalogue Number	Particle Size	Volume
ROXRS025	0.25 - 1.0 mm	0.25 cc
ROXRS05		0.5 cc

#### Storage:

Store at room temperature (15°C/59°F - 30°C/86°F).

RegenerOss® is a registered trademark of Zimmer Biomet or its affiliates.

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