Instruction for use

Zfx™ Cobalt Chromium Eco

Version: 4 / 09.2021
Compiled on: 26.07.2021
Effective date: 01.09.2021
Replaces version: 3 / 10.2018
Page 1 / 2

1. Description

Beryllium and cadmium-free cobalt base alloy CoCrMo6 according to DIN EN ISO 22674 and EN ISO 9693 for applications with thin cross-sections which are exposed to very high forces, e.g. removable partial prostheses, braces, thin and veneered single crowns, fixed full-arch prostheses or bridges with small cross-sections, bars, anchors and implant-supported superstructures

2. Area of application

- × Single crowns
- × Primary and secondary components
- × Bridges in all spans for veneering with dental ceramic
- × Implant construction
- × Fully cast crowns
- × Cone and telescopic technology
- × Fully cast bridges

3. Properties

Type of alloy		4
Density	g/cm ³	8,8
Vickers hardness	HV 5/30	350
CTE value (25-500 °C)	*10 ⁻⁶ K ⁻¹	14,4
Oxidising without vacuum	°C	1.000
Mass content of all parts	%	Co 63, Cr 29, Mo 6, Mn x, Nb x, Si x, Fe x*
Tensile strength	MPa	680
Yield strength 0.2 %	MPa	400
Elongation at break	%	9
E module	MPa	250.000

*x= < 1%

4. Detaching the frameworks

The milled frameworks can be detached from the blank with cross-cut carbide milling cutters suitable for CoCr alloys or with suitable cutting disks.

5. Elaboration/cleaning

The detached frameworks should first be sandblasted with aluminium oxide (approx. 100 µm). Elaboration is carried out with clean carbide milling cutters or diamond grinding wheels suitable for CoCr alloys. Whereby the tool should be moved in only one direction over the surface in order to avoid material overlapping and any resulting blistering during ceramic veneering. In addition the maximum speed of the instruments recommended by the manufacturer must be observed. The surfaces are then sandblasted with pure aluminium oxide (approx. 180 µm) with a pressure of 3 - 4 bar. Next brush off the framework under running water or steam off with superheated steam and degrease with ethanol.

6. Oxidising

Oxide firing is not necessary for CoCrMo6 frameworks. Should oxide firing (rapid heating to 1,000 °C without vacuum with immediate cooling after the temperature is reached) nevertheless be carried out for purposes of visual inspection of the framework conditioning, the oxide must be blasted with disposable blasting abrasive Al₂O₃ (grain size approx. 180 µm) and the surface cleaned again.

7. Firing

All conventional ceramic materials with an appropriate CTE value are suitable for ceramic veneering. The processing instructions recommended by the ceramic manufacturer must be observed, particularly with regard to the cooling conditions after ceramic firing.

Please note:

Due to more favourable cooling conditions in the case of solid pontics it might be preferable to apply a continuous metal collar.

8. Joining

A solder suitable for CoCr alloys is recommended for soldering. The soldering block should be kept as small as possible, and the soldering block must be heated in an oven at 600 °C for approx. 10 minutes before flame soldering is carried out. The soldering surfaces should be brushed with a fluxing agent before heating. Soldering gaps may only be between 0.05 to 0.2 mm. In the case of soldering joints > 0.2 mm there is a risk of the solder shooting through the soldering seam, the joint not being sufficiently filled with solder so that the joint does not have sufficient strength. Allow the soldered part to cool down slowly after soldering. A standard CoCr laser welding wire is recommended in the case of laser welding before ceramic firing.

Attention: Pay attention to the weld seam geometry.

9. Polishing

Remove oxide residues and any fluxing agent residues by mean of blast polishing. Smooth down metal surfaces elaborated with carbide with ceramic-bonded grinding stones and then polish to a semi-matt gloss finish with rubber polishers. Polish with a suitable polishing paste until a high-gloss finish is achieved and then clean with an ultrasonic cleaning device or carefully steam clean.

10. Safety instructions

These medical products may only be used and applied in the dental laboratory by qualified personnel.

Careful: metal vapours and metal dusts are harmful to health if inhaled; for this reason always ensure sufficient extraction and/or wear a suitable protective mask!

11. Explanation of symbols

Manufacturer Information:
Address of legal
manufacturer facility

Distributor Information:
Address of legal
distributor facility

LOT-Number, Symbol followed
by Production Date of the lot.

Medical Device

Reference Number, Symbol followed
by Manufacturer's Item Number

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Our information and recommendations are based on the currently known state-of-the-art in science and technology, and are to be considered correct to the best of our knowledge and experience at this time. This version shall replace all previous versions.

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Wieland Edelmetalle GmbH
Zfx GmbH