# Instruction for use Zfx Titanium G5

Version: 3 / 09.2021  
Compiled at 26.07.2021  
Effective date: 01.09.2021  
Replaces version: 2 / 01.2019  
Page 1 of 2

Ti6Al4V based dental alloy for milling (Grade 5)

Ø 98,5 mm:  
10 mm, 12 mm, 14 mm, 15 mm, 16 mm, 18 mm, 20 mm, 22 mm, 25 mm

## 1. Chemical Composition (Typical Values)

Titan %	Al %	V %	Fe %	O %	C %	N %	H %
Balance	5,5 - 6,5	3,5 - 4,5	0,25	0,13	0,08	0,05	0,012

Determination, extent and validity of the chemical composition according to DIN EN 10 204 – 3.1

## 2. Typical Technical Data (After casting)

Yield strength 0,2 %	816 MPa
Elongation	8 %
Tensile strength	885 MPa
E-module	174 GPa
Density	4,4 g / cm <sup>3</sup>
Corrosion resistance	< 200 µg / cm <sup>2</sup>
Tarnish resistance	ja
Hardness	318 HV 10/30
CTE (25-500°C)	10,3 x 10 <sup>-6</sup> K <sup>-1</sup>
Max. firing temp.	max. 800 °C

## 3. Applied Norm

DIN EN ISO 22674:2016  
ASTM F 136 (z.T.) / DIN EN ISO 5832-3:2012 (z.T.)  
Zfx GmbH is certified according to DIN EN ISO 13485:2016.

## 4. Intended use

Zfx Titanium G5 is a titanium bonding alloy for the production of milled crowns and bridges. Only for professional user! Not intended for the production of implants!

## 5. General guidelines for handling

This instruction for use includes important processing steps and recommendations for Zfx Titanium G5.

## 6. Indication

- Crowns and bridges for the veneering with ceramic and composite with max. 16 units (max. 4 pontics in the anterior region and up to 3 pontics in the posterior region)
- Implant retained superstructures, abutments
- Bars

## 7. Milling

Zfx Titanium G5 is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and CNC milling machine.

## 8. CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6-9 mm2). Sharp edges and undercuts should be avoided.

## 9. Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

## 10. Laser welding

Sufficient and secure connections can only be done by laser welding. The use of Ti-welding wire is an option.

## 11. Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with 110 µm of Aluminium oxide with 2-3 bar and clean with steam cleaner. Consider the firing process with max. 800°C, Otherwise lattice transformation might occur.

## 12. Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

## 13. Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with Titan alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

## 14. Single-use

Used discs should be not processed for further use as a medical device.

## 15. Disposal Instructions

Consult the material safety data sheets or national regulations for disposal. Dispose of Zfx Titanium G5 residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

## 16. Storage conditions

Temperature, humidity or light has no effect on the product properties.

## 17. Quantity

Please consider the label on the package.

## 18. Information on Manufacturer

Manufacturer Information: Address of legal manufacturer facility

Don't use twice

Distributor Information: Address of legal distributor facility

Use by

LOT-Number, Symbol followed by Production Date of the lot.

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

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Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Eisenbacher Dentalwaren ED GmbH