Universal RealGUIDE™ box

IFU RG Box - Rev. 00 – 2023-01

3DIEMME S.r.l., Via Risorgimento 9,
22063 Cantù (CO) – Italy

Contents

| | |
|---|----|
| (EN) English – Instructions for Use | 2 |
| (IT) Italiano – Istruzioni per l'uso | 6 |
| (FR) Français - Instructions d'utilisation | 11 |
| (ES) Español - Instrucciones de uso | 16 |
| (PT) Português - Instruções para o Uso | 21 |
| (DE) Deutsch - Gebrauchsanleitung | 26 |
| (NL) Nederlands - Gebruiksaanwijzing | 31 |
| (BG) Български - Инструкции за употреба | 36 |
| (CS) Čeština - Návod k použití | 41 |
| (DA) Dansk - Brugsanvisning | 45 |
| (HU) Magyar - Használati utasítás | 49 |
| (PL) Polski - Instrukcja obsługi | 53 |
| (RO) Română - Instrucțiuni de utilizare | 58 |
| (SK) Slovenčina - Návod na použitie | 63 |
| (SL) Slovenščina - Navodila za uporabo | 67 |
| (SV) Suomi - Bruksanvisning | 71 |

CE

If any serious incident occurs using the medical device, object of this Instruction for Use, please contact 3DIEMME immediately. Our company will conduct a proper analysis procedure and will report the event to the competent authority of the Member State in which the user and/or patient is established.

1. Product identification

The Universal RealGUIDE™ box is part of the RealGUIDE™ surgical instrumentation for guided surgery. This IFU is applicable to Universal RealGUIDE™ box, REALGUIDE™ Z3D KIT Empty Tray box and Universal RealGUIDE™ MINI box.

The sterilization tray is a hard plastic container, made of 3 parts: a base container attached to the lid with two side hinges and a middle tray that slots in the bottom tray to house the instruments and hold them in place.

Holes on the base container and on the lid allow effective sterilization of the instruments.

An opening handle on the lid allows for easy opening of the device.

Colour codes help the user position the instruments in an orderly manner. Every positioning hole also has the indication of what instrument to position within it.

Materials used allow for the device to be reusable and ideal for withstanding standard steam sterilization cycles.

2. Material

The Universal RealGUIDE™ boxes are made from polyphenylsulfone (PPSU) and silicone. Instrument supports, fixed to the tray are made of stainless steel or Titanium.

3. Intended purpose:

The RealGUIDE™ Sterilization box is a reusable container intended to enclose dental surgery instruments during steam sterilization.

4. Indications for use:

The Sterilization box is indicated to store and organize surgical instruments between and during surgical uses, transport and processing.

It is indicated for use in healthcare facilities by healthcare professionals.

The box is indicated to be used in conjunction with a legally marketed, validated sterilization wrap to maintain the sterility of the devices enclosed.

The surgical instruments must be placed in the silicone holders as indicated on the insert and bottom tray. The surgical instruments requiring connection to the support, must be carefully connected to the sterilization tray as indicated on the insert.

Double pack the box with the instruments in disposable sterilization packaging corresponding to the following requirements (EN ISO 11607):

- Suitable for steam sterilization (temperature resistance up to at least 134 °C (273 °F))
- Sufficient steam permeability
- Sufficient protection of the instruments or sterilization packaging against mechanical damage.

Note: An indicator strip with device ID, LOT. number, the date of sterilization, the expiration date and the operator's name should be affixed to every sterilization packaging.

5. Contra-indications:

None know.

6. Intended users:

The present device must be used only by qualified healthcare professional.

7. Residual risks and undesirable side-effects

No residual risks or undesirable side-effects have been identified.

For any residual risks, contra-indications and any undesirable side-effects related to the surgical instruments, please refer to the specific Instructions for Use, as supplied by the Legal Manufacturer.

8. Cautions/precautions

It is essential to handle the Sterilization Tray according to the basic information documentation for the applicable surgical and prosthetic procedure.

Only use adequately processed devices if indicated for multiple uses.

The Sterilization Tray has been validated to withstand 300 processing cycles provided that the recommended conditions of use are followed. Regardless of the number of times the device has been used, the user must always evaluate its condition after each use.

Corroded or rusted instruments contaminate the water circulation system of the sterilizer with rust particles. During every subsequent sterilization cycle, these rust particles cause rust on instruments that were originally intact. If visible signs of moisture are present (damp spots on sterile packaging, pooled water in the load) at the end of the sterilization cycle, repackage and resterilise using a longer drying time. Make sure that the box and instruments are completely dry before packing for sterilization. Sterilized devices should be used immediately after sterilization. In case of storage, strictly follow the manufacturer's instructions of the sterilization accessories and storage containers.

9. Cleaning

The Sterilization Tray is delivered non-sterile and is not intended to be sterilised before use.

The sterilisation tray must be cleaned prior to every reuse.

For cleaning and disinfection, the tray must be dismantled into its component parts. Do not combine cleaning/disinfectants agents. Always follow the manufacturer's instructions for use for cleaning and disinfection agents.

Always use mild detergents even if cleaning is done using ultrasonic equipment. Dry the boxes with a soft cloth. Do not use abrasive sponges or metal parts that might damage the surface of the box.

10. Handling

Use protective clothing (e.g. protective glasses, facemask, gloves, waterproof gown) while cleaning contaminated instruments.

Take appropriate measures to avoid mechanical damage and contamination, e.g. do not mix heavy devices with delicate ones, safely store and transport the box in its cardboard box.

Pay particular attention to cutting edges, both to avoid injury and to avoid damage to the medical device.

Perform cleaning as soon as possible. When longer delays are expected, immerse the devices in a water bath or disinfection solution to avoid drying of debris.

11. Compatibility information

The Sterilization Tray is intended to hold RealGUIDE™ dental guided surgery instruments during transport, storage and steam sterilization, according to marking and colour coding provided on the tray and base.

Make sure only to use original RealGUIDE™ instruments with the boxes. Storage and organization of non-RealGUIDE™ instruments and components can lead to mechanical and/or instrument failure.

Sterilization must be carried out with disposable packaging, legally marketed and compliant to (EN ISO 11607) requirements.

Cleaning of the tray must be done with a mild detergent.

Sterilization must be carried according to the instructions given by the Instruments' manufacturer.

12. Warnings

Do not use if package is damaged.

Do not use the tray, regardless of the number of cycles of use, if any damage, corrosion, deformation, or cracking is observed, or the surgical workflow becomes illegible. Never use potentially contaminated components.

Do not put the tray on its side or upside down with the lid facing down.

Do not place corroded or rusty instruments in the box for sterilization.

Do not use flash sterilization, hot air sterilization, radiation sterilization, plasma sterilization, formaldehyde, or ethylene oxide sterilization.

Do not sterilize with temperatures higher than 134 °C (273 °F).

Never clean instruments and sterilization boxes with metal brushes or steel wool.

The Sterilization tray does not maintain sterility of the instruments within. Always use appropriate sterilization wraps.

Do not use with other instruments, other than those indicated by the RealGUIDE™ procedure.



13. Storage

The devices must be stored in a dry place at room temperature and be protected from direct sunlight. Improper storage may compromise essential material and design characteristics, leading to device failure.

14. Disposal

Disposal should be handled in an environmentally sustainable manner according to local regulations. Hazardous waste from contaminated devices should be disposed of in appropriate containers which meet specific technical requirements.

15. Symbol glossary

| Symbol | Description | Source |
|---|---|----------------------------|
|  | CE Marking | MDR (EU) 2017/745, Annex V |
| QTY | Contents of the packaging (Quantity of devices) | - |
|  | Reference code (Catalogue number) | ISO 15223-1 |
|  | Batch (LOT) number | ISO 15223-1 |
|  | Legal Manufacturer | ISO 15223-1 |
|  | Keep away from sunlight | ISO 15223-1 |
|  | UDI code | ISO 15223-1 |
|  | Device is Non-Sterile | ISO 15223-1 |
|  | Do not use device if packaging is damaged | ISO 15223-1 |
|  (eIFU) | Read instructions for use (electronic version) | ISO 15223-1 |
|  | Medical Device | ISO 15223-1 |
|  | Authorised Representative in Switzerland | MedDO, Art. 104abis |
|  | Caution: Federal law restricts this device to sale by or on the order of a dentist (USA). | CFR, Title 21, Part 801 |

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