en

## IMPORTANT PRODUCT INFORMATION INSTRUCTIONS FOR USE

# BIOGRAN® 300-355 µm (50-45 mesh)\*

## Resorbable Synthetic Bone Graft

**DESCRIPTION:** Biogran is a resorbable synthetic graft for use in the repair of oral bony defects. It is a bioactive glass (Si, Ca, Na, P) of a narrow size range, which has been shown to accelerate the healing of bone tissue in oral defects.

**ACTIONS:** Biogran is a synthetic bone graft material consisting of 300-355 µm diameter bioactive glass granules. As a result of its optimal size range, Biogran is transformed into hollow calcium phosphate bone growth chambers. Each bone growth chamber within a bony defect provides a protective environment in which cells differentiate into osteoblasts, the producers of new bone. Osteogenesis occurs at multiple sites, rapidly filling the osseous defect with new bone that continuously remodels in the normal physiological manner.

**INDICATIONS:** Clinical Studies have demonstrated the effectiveness of Biogran for treating periodontal defects, extractions sites and reconstructions and augmentation of the alveolar ridge.

**CONTRAINDICATIONS:** Bone grafting should not be considered for patients where general oral surgery is not advisable.

**PRECAUTIONS:** Underlying oral pathological conditions, e.g. infections, should be controlled or eliminated prior to implementation of Biogran. PRIMARY CLOSURE OF THE DEFECTS MUST BE ACHIEVED. Care should be exercised to avoid salivary contamination.

Reuse of BIOMET 3i products that are labeled for single-use may result in product contamination, patient infection and/or failure of the device to perform as intended.

**ADVERSE REACTIONS:** No adverse reactions have been reported.

**HOW SUPPLIED – AVAILABLE IN DIFFERENT SIZE DOSES:** **Cups** – Biogran granules are conveniently packaged in individually sealed dappen dishes. **Syringes** – Glass syringes are packaged in individually sealed trays.

**INSTRUCTIONS FOR USE: Cups** – The seal on the cup should be carefully peeled back to avoid spillage. Biogran can be mixed with the patient's blood or sterile saline within the disposable dappen dish cup and subsequently delivered to the defect site with a dental spatula. **Syringes** – The syringe package offers a choice of delivery methods. Biogran granules can be delivered directly to the prepared oral bone defect site from the syringe prior to soft tissue closure. Mix the Biogran with the patient's blood or sterile saline within the disposable dappen dish. Pull back on the syringe plunger to create a space for the mixed Biogran. Pack the syringe cylinder by tamping the open end of the syringe into the mixed Biogran in the dappen dish. Continue packing until the volume of Biogran you wish to dispense is in the syringe cylinder. Carry the syringe to the defect site with the syringe tip pointed upward. Deliver the Biogran to the surgical site by pushing the plunger forward. Biogran can also be pre-mixed with the patient's blood or sterile saline within the disposable dappen dish cup and subsequently delivered to the defect site with a dental spatula.

**GAMMA STERILIZED:** Syringe, cup interior and contents are sterile unless lid is opened or damaged. Discard unused portion immediately after use. For single use only.

**Caution:** U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

**LOT** Batch code

**REF** Catalogue number

 Caution, consult accompanying documents

 By prescription only

 Do not re-sterilize

 Do not re-use

**STERILE R** Sterilized using gamma irradiation

 Use by date

 Do not use if package is damaged

 Consult Instructions for Use

**BIOMET 3i**

BIOMET 3i  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454  
Outside The U.S.: +1-561-776-6700  
Fax: +1-561-776-1272  
www.biomet3i.com

**EC REP**

BIOMET 3i Dental Iberica, S.L.  
WTC Alameda Park, Ed. 1, Planta 1ª  
Pl. de la Pau s/n  
08940 – Cornellà de Llobregat  
(Barcelona) Spain  
Phone: +34 934 705 500  
Fax: +34 933 717 849

Biogran is a registered trademark of BIOMET 3i LLC. BIOMET is a registered trademark and BIOMET 3i and design are trademarks of BIOMET, Inc. ©2008 BIOMET 3i LLC. All rights reserved. \*ISO3310, ASTM E11

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