en

**IMPORTANT PRODUCT INFORMATION  
INSTRUCTIONS FOR USE**

**BIOGRAN®** 300-355 µm (50-45 mesh)\*

**Resorbable Synthetic Bone Graft**

**DESCRIPTION:** Biogran is a resorbable synthetic graft for use in the repair of oral bony defects. It is a bioactive glass (Si, Ca, Na, P) of a narrow size range, which has been shown to accelerate the healing of bone tissue in oral defects.

**ACTIONS:** Biogran is a synthetic bone graft material consisting of 300-355 µm diameter bioactive glass granules. As a result of its optimal size range, Biogran is transformed into hollow calcium phosphate bone growth chambers. Each bone growth chamber within a bony defect provides a protective environment in which cells differentiate into osteoblasts, the producers of new bone. Osteogenesis occurs at multiple sites, rapidly filling the osseous defect with new bone that continuously remodels in the normal physiological manner.

**INDICATIONS:** Clinical Studies have demonstrated the effectiveness of Biogran for treating periodontal defects, extractions sites and reconstructions and augmentation of the alveolar ridge.

**CONTRAINDICATIONS:** Bone grafting should not be considered for patients where general oral surgery is not advisable.

**PRECAUTIONS:** Underlying oral pathological conditions, e.g. infections, should be controlled or eliminated prior to implementation of Biogran. PRIMARY CLOSURE OF THE DEFECTS MUST BE ACHIEVED. Care should be exercised to avoid salivary contamination.

**ADVERSE REACTIONS:** No adverse reactions have been reported.

**HOW SUPPLIED – AVAILABLE IN DIFFERENT SIZE DOSES:** **Cups** – Biogran granules are conveniently packaged in individually sealed dappen dishes. **Syringes** – Glass syringes are packaged in individually sealed trays.

**INSTRUCTIONS FOR USE:** **Cups** – The seal on the cap should be carefully peeled back to avoid spillage. Biogran can be mixed with the patient's blood or sterile saline within the disposable dappen dish cup and subsequently delivered to the defect site with a dental spatula. **Syringes** – The syringe package offers a choice of delivery methods. Biogran granules can be delivered directly to the prepared oral bone defect site from the syringe prior to soft tissue closure. Mix the Biogran with the patient's blood or sterile saline within the disposable dappen dish. Pull back on the syringe plunger to create a space for the mixed Biogran. Pack the syringe cylinder by tamping the open end of the syringe into the mixed Biogran in the dappen dish. Continue packing until the volume of Biogran you wish to dispense is in the syringe cylinder. Carry the syringe to the defect site with the syringe tip pointed upward. Deliver the Biogran to the surgical site by pushing the plunger forward. Biogran can also be pre-mixed with the patient's blood or sterile saline within the disposable dappen dish cup and subsequently delivered to the defect site with a dental spatula.

**GAMMA STERILIZED:** Syringe, cup interior and contents are sterile unless lid is opened or damaged. Discard unused portion immediately after use. For single use only.

**Caution:** U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

<b>LOT</b> Batch code	<b>REF</b> Catalogue number
Caution, consult accompanying documents	By prescription only
Do not re-sterilize	Do not re-use
<b>STERILE R</b> Sterilized using gamma irradiation	Use by date
Do not use if package is damaged	Consult Instructions for Use

**BIOMET 3i**  
 BIOMET 3i  
 4555 Riverside Drive  
 Palm Beach Gardens, FL 33410  
 1-800-342-5454  
 Outside The U.S.: +1-561-776-6700  
 Fax: +1-561-776-1272  
 www.biomet3i.com

**EC REP**  
 BIOMET 3i Dental Iberica, S.L.  
 WTC Alameda Park, Ed. 1, Planta 1ª  
 Pl. de la Pau s/n  
 08940 – Cornellà de Llobregat  
 (Barcelona) Spain  
 Phone: +34 934 705 500  
 Fax: +34 933 717 849

**CE**  
**0086**

P-BGS  
 Rev. D

BIOMET 3i and Biogran are registered trademarks of BIOMET 3i, ISO3310, ASTM E11

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