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### **CopiOs® Pericardium Membrane**

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Please read carefully!

### CopiOs® Pericardium Membrane

Collagen membrane for guided bone and soft tissue regeneration

**STERILE – single use for one patient only**

#### Product name

CopiOs Pericardium Membrane

#### Description

The CopiOs Pericardium Membrane is a collagen membrane from solvent-preserved bovine pericardium sterilised by gamma irradiation.

The tissue is selected according to the strictest specifications. CopiOs Pericardium Membrane is preserved using the company's proprietary multiple step process (Tutoplast® process), which consists of a thorough purification of the tissue and gentle solvent dehydration; the process is subject to continuous quality controls. Preservation also includes sterilization of the finished product with low-dose gamma irradiation. The Tutoplast process maintains the natural collagen structure and thus the high biomechanical stability of the tissue.

#### Mechanism of action

The CopiOs Pericardium Membrane acts as a template for the local regeneration of vascularised connective tissue. When applied between bone graft material and soft tissues, the CopiOs Pericardium Membrane takes on the functions of a barrier membrane.

#### Indications

- Guided tissue regeneration (GTR) in periodontology
- Covering and protecting bone graft material, e.g. in guided bone regeneration (GBR)

#### Contraindications

Known hypersensitivity to bovine collagen.

Disorders or diseases that carry an unacceptable increase in postoperative risk.

#### Warnings

The device should be implanted adjacent to well-vascularized tissue for best infiltration of cells, revascularization and incorporation.

CopiOs Pericardium Membrane remains sterile, as long as the sterile packaging is not opened or damaged. Sterility is no longer ensured once the package has been opened. If the product is resterilised, the stability/quality of the material can no longer be guaranteed. Re-use is therefore not permitted as it may harm the patient's health.

Integration and remodeling can be affected by comorbidities or conditions that lead to a reduced vascularisation, such as diabetes or smoking.

#### Precautions

Use caution when implanting in areas that are infected or contaminated. The most recent clinical reports on individual cases allow the assumption that the use of xenogeneic membranes such as CopiOs Pericardium Membrane is safe and successful under such conditions. However, there are not sufficient clinical data to support the routine use of CopiOs Pericardium Membrane in infected, contaminated or possibly infected or contaminated areas. A careful risk assessment must therefore be made for each individual patient.

Inspect the device, including all packaging and labeling materials carefully:

- Do not use past expiration date.
- Do not use if the device or packaging is damaged.
- Do not use if there are discrepancies in label information.
- Do not use if the device was rehydrated in any solution above 37°C.

#### Handling instructions

Please ensure that the general rules for handling medical devices under aseptic conditions are observed when using CopiOs Pericardium Membrane.

#### Special precautions for use

The Tutoplast process is highly effective against all pathogenic species. However, as with all biological products, the transmission of infectious diseases can be ruled out only to the extent of current knowledge.

## Preparation

- Select the right CopiOs Pericardium Membrane to provide adequate cover. CopiOs Pericardium Membrane does not shrink significantly during the integration process.
- Remove the pouch from the carton.
- Peel open the outer pouch and remove the inner pouch using aseptic technique.
- The inner pouch is sterile and can be placed directly into the sterile field.
- Open the inner pouch and aseptically remove the device with sterile gloves or forceps. Place the tissue on a sterile surface.
- The CopiOs Pericardium Membrane should be cut while still dry and placed over the defect.
- The device must be moist at all times. Rehydrate CopiOs Pericardium Membrane with sterile normal saline or Ringer solution (max. temperature 37 °C) or the patient's blood until the membrane is soft and flexible and thus well adapted to the surface. Keep the device in room temperature sterile saline or Ringer solution until ready for implantation.
- Place the CopiOs Pericardium Membrane where it is required and smooth it out from the centre towards the edges so that it lies firmly on the underlying tissue base. The selection of the size should ensure a tension-free placement with sufficient overlap to healthy tissue. The membrane should be adapted to the size of the defect and positioned securely to prevent displacement and to encourage its incorporation. The CopiOs Pericardium Membrane should be implanted in such a way that the free edges do not protrude. Every effort should be made to ensure a tension-free hermetic flap closure.

## Implantation

- Using standard techniques, prepare the surgical site.
- Transfer the device to the surgical site.
- The device may be trimmed prior to fixation with blunt scissors.
- Fix to the patient's tissues with appropriate material with a round-bodied atraumatic needle. If using sutures, they should be placed at least 5 mm from the margin of the collagen membrane.
- Unused CopiOs Pericardium Membrane must be disposed of correctly.
- If dehiscence occurs, the exposed CopiOs Pericardium Membrane does not necessarily have to be explanted. The wound should be disinfected daily (e.g. chlorhexidine) until the membrane is covered with mucosa.
- Patients should be instructed to keep a good dental hygiene, as this is a basic requirement for a successful augmentation.
- Absorption and tissue regeneration processes begin after one or two days and may last for weeks or even months.

## Interactions

None known

Magnetic Resonance (MR) environment: As the device is composed of natural collagen, it is safe for magnetic resonance environment. No MR artefacts are expected due to the composition of the material.

## Adverse reactions

Local reactions may rarely occur in patients who are not known to be hypersensitive to foreign collagen; however, these reactions resolve rapidly. Patients should be told to inform the doctor about any possible side effects. All adverse reactions must be reported to Tutogen Medical GmbH.

## Shelf life

Do not use the CopiOs Pericardium Membrane after the expiry date (  ) printed on the inner pouch and cardboard box.

## Storage instructions

Store this medical device at room temperature between 15 °C and 30 °C, in a clean, dry place; protected from direct sunlight.

## Date of last revision

"Draft"

CopiOs® is a registered trademark of Zimmer Biomet or its affiliates.  
Tutoplast® is a registered trademark of Tutogen Medical GmbH.

**Manufacturer:**

Tutogen Medical GmbH  
 Industriestrasse 6  
 91077 Neunkirchen am Brand  
 Germany  
 Tel.: +49 (0) 9134 99 88-0  
 Fax: +49 (0) 9134 99 88-99  
 tutogen@rtix.com

**Distributed by:**

ZIMMER BIOMET

BIOMET 3i Dental Ibérica, S.L.  
 WTC Almeda Park, Ed. 4, Pl. 2  
 C/Tirso de Molina, 40  
 08940, Cornellà de Llobregat, Spain  
 Tel: +34 93 470 59 50  
 Fax: +34 93 372 11 25  
 zimmerbiometdental.com

 Do not reuse	 Manufacturer	 Keep away from direct sunlight
 Do not resterilise	 Sterilised by irradiation	 Keep in a dry place
 Batch number	 Use by	 Storage temperature
 Order number	 Read the instructions for use	 Do not use if the packaging is damaged

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Tutoplast® es una marca comercial registrada de Tutogen Medical GmbH.