
MS303-2 Rev 2017

**Labeling Symbols:**  
 Symbols may be used on some international package labeling for easy identification.

 Caution, consult accompanying documents

 Catalog Number

 Expiration Date

 USA law restricts this device to sale by or on the order of physician

 Do not reuse after opening

 Manufacturer

 Lot number

 Authorized European Representative

 Method of sterilization – gamma

 Do not use after the product sterilization barrier or its packaging is compromised.

 Temperature Limitation

**Zimmer®**  
**Socket Repair Membrane**

**ENGLISH**

#### Zimmer® Socket Repair Membrane Instructions For Use

Distributed by:  
 Zimmer Biomet Dental  
 4555 Riverside Drive  
 Palm Beach Gardens  
 FL 33410 USA  
 +1.561.776.6700  
 zimmerbiomet.com

#### Intended Use:

*Zimmer®* Socket Repair Membrane is a bioresorbable, implantable collagen material that is intended for use in dental surgery procedures as a material for placement in the area of dental implant, bone defect or ridge reconstruction to aid in wound healing post dental surgery.

#### Description:

*Zimmer®* Socket Repair Membrane is a white, nonfriable membrane matrix engineered from highly purified type I collagen fibers derived from bovine Achilles tendon. *Zimmer®* Socket Repair Membrane is resorbable which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane.

*Zimmer®* Socket Repair Membrane has a morphology of dense oriented fibers for mechanical strength. Macromolecular permeation studies have shown that the membrane is permeable to macromolecules. Its porosity is such that it effectively retards epithelial down growth and prevents gingival connective cell migration into the wound site. The semi-permeability properties of the membrane permit the exchange of essential nutrients for wound healing.

*Zimmer®* Socket Repair Membrane is supplied sterile, non-pyrogenic and for single use only.

**Administration:**  
*Zimmer®* Socket Repair Membrane is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The periodontal defect or bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). The clinician should perform thorough debridement and good planning of the defect. Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect. As much tissue as possible should be preserved to allow for primary closure of the wound and correct positioning of the flaps.

*Zimmer®* Socket Repair Membrane can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to the final placement.

*Zimmer®* Socket Repair Membrane can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

*Zimmer®* Socket Repair Membrane should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. An attempt should be made to suture the mucoperiosteal flap over the collagen membrane to avoid accelerated resorption due to membrane exposure.

#### Post-operative Procedures:

*Zimmer®* Socket Repair Membrane is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

*Zimmer®* Socket Repair Membrane should be completely resorbed 26 to 38 weeks following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

#### Contraindications:

*Zimmer®* Socket Repair Membrane is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity
- known allergy to collagen of animal origin or other bovine-derived products
- clinically significant renal, hepatic, cardiac, endocrine, hematological, autoimmune or systemic disease, which in the physician’s judgment, will prevent safe implantation or likely healing.

#### Warning:

Clinicians should use care in screening their patients for any known allergies to collagen or bovine-derived products. Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to *Zimmer®* Socket Repair Membrane.

#### Precautions:

As with all surgical procedures, caution should be exercised when treating medically compromised patients. In this study, the following steps were evaluated for their ability to inactivate or reduce the following viral strains: Bovine Viral Diarrhea (enveloped virus) and Porcine Parvoviridae (non-enveloped virus). The study results showed that each of the manufacturing steps evaluated, including the sodium hydroxide treatment, is effective in inactivating these viruses.

*Zimmer®* Socket Repair Membrane cannot be resterilized. Open, unused *Zimmer®* Socket Repair Membrane must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

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**Safety:**  
 The product is manufactured from bovine Achilles tendon, which is classified as tissues with no detected infectivity for Bovine Spongiform Encephalopathy, BSE (World Health Organization Guidelines). The bovine tendon is known to be one of the richest sources of type I collagen that is commercially available.

The manufacturing process for the product meets European Standards and International Standards for animal tissue sourcing, handling and inactivation of Spongiform Encephalopathy (SE) pathogens. This process involves a treatment with sodium hydroxide that is a recognized method of inactivation of SE pathogens.

A viral inactivation study for the product’s manufacturing process was conducted by an independent laboratory. In this study, the following steps were evaluated for their ability to inactivate the following viral strains: Bovine Viral Diarrhea (enveloped virus) and Porcine Parvoviridae (non-enveloped virus). The study results showed that each of the manufacturing steps evaluated, including the sodium hydroxide treatment, is effective in inactivating these viruses.

#### Storage:

The product should be stored at room temperature. Avoid excessive heat and humidity.

#### How Supplied:

One (1) membrane per package, size 10mm x 20mm

#### Caution:

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

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©2017 Zimmer Biomet. Todos los derechos reservados. Zimmer es una marca registrada de Zimmer Biomet o sus filiales.
**Attenzione:** secondo la legge federale (U.S.A.) questo prodotto deve essere venduto solo a medici o dentisti.